






# Greenwatch v2 GUI Flow and Design Document

Revision 1

Prepared by

Sharome Burton  
Delton Hughes  
Victor Marchesi  
Calvin Walmer

March 26, 2024



To the Room page

Create User

Edit Room Popup

Create User Popup

Edit User Popup

From Create Room Button

From Create Room Button

↓

Create Room

Room Name

Required fields cannot be empty.

Cancel Create

From Edit Room Button

From Edit Room Button

↓

Edit Room

Editing: Room

Room Name

Room

Cancel Save

From Delete User Button

From Delete User Button

↓

Create User

First Name

Last Name

Username

Password

Email

Admin Status

Select Status

Required and:  
- contains one uppercase letter  
- contains one lowercase letter  
- contains one number  
- contains one special character  
- min 8 to max 100000000

Cancel Create

From Edit User Button

From Edit User Button

↓

Edit User

Editing: user 1

First Name

Last Name

Username

Password

Email

Admin Status

Select Status

Cancel Save

From Roompage

**Room Page**

**Bolin [10.20.52.40]**

Start Date: 03/25/2024 End Date: 03/25/2024

Update Date Range Export Data

Chart View: All

Temperature Humidity Pressure Light

Experiments

	Measurement	Low	High
<input checked="" type="checkbox"/>	Temperature	1	50
<input checked="" type="checkbox"/>	Humidity	1	90

**Start Date**

Start Date: 03/25/2024 End Date: 03/25/2024

Update Date Range Export Data

Chart View: All

Temperature Humidity Pressure Light

**End Date**

Start Date: 03/25/2024 End Date: 03/25/2024

Update Date Range Export Data

Chart View: All

Temperature Humidity Pressure Light

**Chart Selection**

Start Date: 03/25/2024 End Date: 03/25/2024

Update Date Range Export Data

Chart View: All

Temperature Humidity Pressure Light

**Experiment**

New Experiment

Experiment Name

Start Date: mm/dd/yyyy

End Date: mm/dd/yyyy

Temperature High

Temperature Low

Humidity High

Humidity Low

Input fields cannot be empty.

Cancel Add Experiment